

Automatic Blood Bank Tube Sealer

Model: GR2002

The tube sealer is newly developed for sealing the blood bag. It is widely used in the blood bank center with more reliable efforts and a more attractive appearance.



Technical data:

Sealing Voltage	adjust automatically
Control	microprocessor control
Sealing Time	1 ~ 2s
Size of Tube	Φ4mm to Φ6mm thickness ≤ 1mm
Power Supply	AC220V 50/60Hz 2A
NW	7 kg
GW	8 kg
Dimension	280 x 115 x 220 mm
Packing Size	450 x 170 x 280 mm